

## ASSEMBLEE GOO Firm and Light Tack

World's Finest Lubricants!

Dr. Tranny's

Assemblee

GOO

Transmission Assembly Lube

Net Wit To be (Ind/454 grams)

#19250



#19260

Product	Dr. Tranny #19250	Dr. Tranny #19260
Melting Temperature	130°F	116°F
Adhesive Strength	1.1 lbs.	1.1 lbs.
Color	Green	Blue
Rust Inhibitor Added	Yes	No

**DESCRIPTION:** Dr. Tranny™ Assemblee Goo™ is formulated for use as an assembly aid and lubricant for automatic and manual transmissions. It provides lubrications and helps facilitate the assembly process, by holding components such as needle bearings and hydraulic seals in place. Simply apply to O-rings, seals, bearings, gaskets, sealing rings, bushings, thrust washers and much more.

Dr. Tranny Assemblee Goo remains smooth and stable. It will not melt at shop or work bench temperatures and provides either light or firm tack for your assembly requirements. Yet it dissolves and becomes soluble in transmission fluids at operating temperatures without clogging filters or solenoid valves. Dr. Tranny Assemblee Goo eliminates misdiagnosing leaks and is NOT harmful with prolonged skin contact.

## Now available in two tack types for the professional technicians preference:

GREEN: DR.TRANNY's original firm tack formula (GREAT FOR SUMMER) BLUE: DR.TRANNY's light tack formula (GREAT FOR WINTER)

## **BENEFITS:**

- Provides your choice of tack strength, firmer (GREEN) great for summer months and lighter tack (BLUE) great for winter months.
- Facilitates and provides lubrication during the transmission assembly process
- Dyed green or blue to prevent misdiagnosing leaks
- Will not melt at shop temperatures
- Not harmful with prolonged skin contact
- Will not clog filters
- Dissolves quickly at operating temperatures
- Compatible with all transmission fluids

**DIRECTIONS:** Apply a dab of Dr. Tranny Assemblee Goo for use as an effective aid in holding components in place during the assembly process. It is an assembly lubricant for all automatic and manual transmissions.

## **AVAILABILITY:**

Stock No.	Unit Size	Case Weight/QNTY
19250 (Green)	16 oz.	13.58 lbs. (12/c)
19260 (Blue)	16 oz.	13.36 lbs. (12/c)





7930 Occidental South Seattle, WA 98108 PO Box 24743 Seattle, WA 98124 206-762-5343 800-333-LUBE (5823) Fax: 206-762-7989 www.Lubegard.com